

Xdigi-USB

IOPA X-Ray Digitizer



Product Specification
 155 mm Height
 85 mm Width
 50mm Depth
 300 gms Weight

Now with updated Interface Software



Product Details :

- * Simplest way to digitize your IOPA X-Ray.
- * Simple 5 Step Process :
 Shoot X-Ray, Develop & fix, put it in Xdigi-USB device and just grab the X-Ray.
- * Convert your physical X-Ray into digital format.
- * Must be tool for better Patient Explanation.
- * USB Interface to connect to Laptop / Desktop.
- * Image Management & Interface software included.
- * More than 15 image enhancement tools included.
- * Brightness/Contrast & Length Measurement at its best.
- * Image Transfer & Exchange Facility.
- * Can be integrated with Patient Management Software.
- * Can be used to convert all your old IOPA X-Ray for storage.
- * Small in size and can be transferred very easily.
- * Compatible with all versions of Microsoft Windows.
- * Interface software has features like, Brightness/Contrast, Zooming, Spot magnification, Canal Measurements, uneven area markings, Flip & Mirror, Crop, Auto Contrast, Emboss and more...
- * Regular Software up-gradation from net with added features.


Xdigi-USB X-Ray Digitizing Process



Shoot IOPA X-Ray on Standard Plate.



Develop & Fix your X-Ray.




Put Developed X-Ray on Glass Plate.



Put the glass plate in Xdigi-USB Device & Grab the X-Ray.



Your digital version of IOPA X-Ray is ready & enjoy Image Editing functionality on it.



Special Price
₹ 8,500/-



Designed & Developed By :
 Jolly Computers,
 Regd. Office : 309, Dalal Chambers,
 Nr. Telephone Exchange, Nanpura,
 Surat - 395 001. Gujarat INDIA.
 Tele / Fax : (0261) 2476427.
 Mobile : 098241 34266. RIM : 093746 65580
 Email : jollycom@yahoo.com / info@jcindia.com
 URL : www.jcindia.com